


TOUCHLESS AUTOMATIC IR THERMOMETER



DhonaadhiTM
Hitec Innovations



Call : +91 98841 66566

Mail : Info@dhonaadhi.in

Web : www.dhonaadhi.in

**Thank you for purchasing our product , Please Read the USER MANUAL
carefully before Operating.**



Accessories



Thermometer



1800MM Cable



18650 Battery

Optional Accessories




Adapter




Tripod



FEATURES



Measuring distance





Accessories


1. Please connect the power supply to activate the device when the battery is installed for the first time.
2. Battery charging state. a GREEN light means 'fully charged'. RED LIGHT Means 'in-charging'. FLASHING LIGHT Means ' LOW BATTERY '
3. There are three power sources available for the device.
 - a. 18650 Li-ion rechargeable battery.
 - b. Power supply through the USB cable.
 - c. Any Power bank as an external battery backup by connecting it to the USB port.

STATUS DESCRIPTION

1. Standby : Long press the power button for 3 seconds until you hear the "beep" sound. The red dot in the bottom of the display turns on indicating Standy mode.



2. Normal temperture : Device beeps once with green flash indicating normal temperature.
3. Elevedaed temperature : Device beeps thrice with red flash indicating elevated temperature.
4. Modes :
 - a. Bod mode (default)
 - b. Ultra-low ("Lo") temperature $< 35^{\circ}\text{C}$
 - c. Ultra-high ("Hi") temperature $> 42^{\circ}\text{C}$






OPERATING INSTRUCTION

There are three buttons in the centre of the device. “+”, “power / setting”, “-” (from left to right).




1. Switch ON/OFF the Device

- ON : Long press “power/setting” button for 3 seconds, all lights in the display will turn on.
- OFF : Long press “power/setting” button for 3 seconds. all the lights in the display will turn off.



2. MODE Selection

- Long press “+” button for 3 seconds
- SUR : Mode for measuring object temperature. Temperature range is 0°C - 100°C
- COU : Mode for counting. Counting capacity . 0-9999
- BOD : Mode for measuring body temperature. Temperature range is 34°C - 43°C



3. Temperature unit switching between “°C / °F ”.

- Press both “+” and “-” button together until the display shows “F1”.
- Select “+” for “C and ”-” for °F.



4. Reset Alarm Temperature

- Press both “+” and “-” button together until the display shows “F1”




- b. Then press the power/setting button and the display will show F2.
- c. Press “+” or “-” to reset alarm temperature.
- d. Press Power button until you hear a long beep.



5. Temperature Devision Setting

- a. Press both “+” and “-” button together until the display shows “F1”.
- b. Then press the power/setting button twice and the display will show “F3”.
- c. Press “+” or “-” to increase / decrease temperature (the difference value is 0.1°C). If the season or the environment changes, the equipment should be checked and adjusted.



6. Mute (Disable) Alarm

- a. Press both “+” and “-” button together until the display shows “F1”.
- b. Then press the power / setting button thrice and the display will show “F4”.
- c. Then Press “-” to disable the alarm and press “+” to enable the alarm



**SPECIFICATIONS**

Accuracy	$\pm 0.2^{\circ}\text{C}$ (34-45$^{\circ}\text{C}$, Place it in the operating environment for 30 minutes before use)
Abnormal Automatic Alarm	Flashing + “Di Di” sound
Automatic Measurement	Measuring distance 5cm- 10cm
Screen	Digital display
Charging Method	USB charging or battery (18650 Li-ion)
Installation Method	Nail hook, double - sided adhesive sticking, Bracket fixing, or Screw to Tripod
Environment Temperature	10$^{\circ}\text{C}$ - 40$^{\circ}\text{C}$ (Recommanded 15$^{\circ}\text{C}$ - 35$^{\circ}\text{C}$)
Infrared Measuring Range	0$^{\circ}\text{C}$ - 50$^{\circ}\text{C}$
Response Time	0.5 Seconds
Input	5V DC
Weight	380 Gram
Diamensions	170 x 115 x 140 mm
Stand by	72 hours
Counting Capacity	9999



Note : After entering the settings menu, if you need to modify a setting , then long-press the menu to select. Save the new settings, you will hear a long beep which indicates that the new settings are saved. If the device Shuts down and exits automatically it indicated that the new settings are not saved.

Incase any queries or issues , Please feel free Mail Us Info@dhonaadhi.in